

- [1] Mackenzie I, Xu J, Cusick C, Widwinter-Morton H, Meacher H, Mollison J, et al. Acupuncture for pain relief during induced labor in nulliparae: a randomized controlled study. *Obstet Anesth Dig* 2012; 32: 52-3.
- [2] Pan P, Bogard T, Owen M. Incidence and characteristics of failures in obstetric neuraxial analgesia and anesthesia: a retrospective analysis of 19,259 deliveries. *Int J Obstet Anesth* 2004;13:227-33.
- [3] Hermanides J, Hollmann M, Stevens M, Lirk P. Failed epidural: causes and management. *Br J Anaesth* 2012; 109:144-54.
- [4] Motamed C, Farhat F, Rémérand F, Stéphanazzi J, Laplanche A, Jayr C. An analysis of postoperative epidural analgesia failure by computed tomography epidurography. *Anesth Analg* 2006;103:1026-32.
- [5] Tielens LK, Bruhn J, Vogt M, van Geffen GJ, Scheffer GJ. The Episure Autodetect syringe, a loss-of-resistance technique for locating the epidural space, used in pediatric patients. *Paediatr Anaesth* 2013; 23:747-50.
- [6] Levin J, Wetzel R, Smuck M. The importance of image guidance during epidural injections: rates of incorrect needle placement during non-image guided epidural injections. *J Spine* 2012;1:1-3.
- [7] Eappen S, Blinn A, Segal S. Incidence of epidural catheter replacement in parturients: a retrospective chart review. *Int J Obstet Anesth* 1998;7:220-5.
- [8] Xiaofeng Liu, E-er-dun Wang, Qing Yan, Kezhong Li. Clinical application of a novel developed pressure bladder indicator in lumbar epidural puncture. *Journal of Clinical Anesthesia*, 2015;27:543-5.
- [9] Jize Feng, Jincai Zou, Jianhua Wong. The application value of LOR Indicator Syringe in epidural puncture. *China Journal of Pharmaceutical Economics*, 2015;6:125-6.
- [10] Junmei Hou, Wenping Jia, Feifei Zhao. The effectiveness comparison of LOR Indicator Syringe with positive pressure syringe in epidural puncture. *Hebei Medical Journal*, 2016;08:1234-6.

TUORen
To Be Better



Tel/Fax: +86-373- 8747537

Email: info@bjtuoren.com; info@etuoren.com

Add: Tuoren Industrial Zone, Changyuan County, Henan Province, China



Epidural Puncture **LOR Indicator Syringe**



TUORen

To Be Better



P Clinical problems

According to Mackenzie et al [1], epidural anesthesia is one of the most difficult procedures to perform in anesthesiology. High failure rate occurs in epidural anesthesia.

- According to the study of Pan et al, a failure rate of 14% was found in epidural anesthesia for laboring women [2].
- Hermanides et al [3] also investigated a number of studies on epidural anesthesia, and the overall failure rate was about 30%. Among the failed epidurals, incorrect catheter placement related to identification of epidural space was responsible for half of the failures, Motamed said [4].

Puncture Maneuver

Correct placement obviously requires correct identification of the epidural space. While, the most common method used for epidural space identification is sudden loss-of-resistance (LOR) when the tip of the needle passes the ligamentum flavum. (Tielens et al [5])

The traditional methods are overdependent on experiences and skills of anesthetists. (Levin J et al [6])



For inexperienced practitioners

- It's very difficult to learn and judge the epidural space by feeling. To the physicians, it is hard to see whether their internships puncture to the right position or not, increasing the difficulty of clinical teaching.
- The loss-of-resistance is a subjective feeling, higher failure rates occur with inexperienced practitioners. (Eappen et al [7])

For experienced anesthesiologists

Even with rich experience, they cannot have a visible judgment and repeatedly pushing action and feeling make the puncture process more complex.



No. 401 Hospital of PLA (Chinese People's Liberation Army) has developed a pressure bladder indicator, which can provide a visible “brake signal” to indicate correctly that the needle has reached the epidural space. Tuoren Medical transforms the technique into a product named LOR Indicator Syringe.

- **Intended Use**

LOR Indicator Syringe is used for epidural anesthesia, combined spinal-epidural anesthesia to indicate that the epidural needle reaches the epidural space.


- **Principle**

The pressure in balloon is higher than that of epidural space, so the balloon gets deflated immediately when the needle reaches the epidural space.



Notice

Sometimes, the balloon gets deflated slowly because not the whole needle tip enters the epidural space. In this case, advancing the needle 1 to 2mm is suggested.



- Confirm the LOR Indicator Syringe is intact before operation; Fill the syringe with about 5 ml normal saline; puncture the epidural needle into dense tissue.



- Remove the stylet, and connect the prepared syringe with epidural needle tightly.



- Push the plunger rod to the mark of 2ml to make the balloon inflated; Next, push the needle forward slowly, and avoid moving backwards.



- When the balloon gets deflated promptly, it means the needle reaches the epidural space successfully.

Accuracy and reliability

The combination of visible indication and experience improves the accuracy and reliability to identify the epidural space.

Safety

The pressure in inflated balloon is 9 ± 2 KPA, which makes the puncture process safer.

Simple operation

Reduce the actions of pushing the syringe repeatedly. Anesthesiologists can control the needle with both hands, improving the stability and effectiveness.

Low resistance

Have the same low resistance function as the traditional LOR syringe.

Visibility

Compared with the blinded traditional syringe, LOR Indicator Syringe makes the puncture process visualized.

Clinical teaching

The visibility can allow the teacher and student to judge the depth of the needle tip intuitively, which facilitates clinical teaching and improves the success rate of the epidural administration.



Study & Research

The balloon pressure is 9 ± 2 KPA.
Why?


- Too low pressure in the balloon might not be able to deflate the balloon when the epidural needle reaches the epidural space. (false-negative indication)
- Extremely high working pressure would make balloon shriveled before the needle reaches the epidural space. (false-positive indication)



Test Data

Two groups were tested (30pcs LOR Indicator Syringes for each group). In the first group, syringes are filled with 3ml normal saline, the second 5ml normal saline.


Hydrapress Measurement of Pressure-indicated Balloon



The novel developed pressure bladder indicator was a reliable and useful technique to conduct successful lumbar epidural puncture [8].

High Success Rate

- Yan et al [8] reported that 70 patients who were undergoing lumbar epidural anesthesia or combined spinal-epidural anesthesia were enrolled to use LOR Indicator Syringe. The success ratio was 100%, indicating that the sensitivity of the bladder indicator in confirming the correct location of the epidural needle is high.
- According to the study of Jize et al, 368 patients accepted epidural anesthesia with LOR Indicator Syringe. The success rate of epidural catheter placement was 100% [9].
- In the essay of Hou et al [10], the success rate and effectiveness of puncture and epidural catheter placement to use LOR Indicator Syringe is higher than using traditional syringe.



Name	Type	Quantity per box	Quantity per carton
LOR Indicator Syringe	5ml	20pcs	400pcs

KITS WITH LOR INDICATOR SYRINGE



	Combined Spinal & Epidural Kit				Continuous Epidural Kit			
Pro. Cod	SE11625C	SE11625R	SE11625T	SE11827C	CE11690C	CE11690R	CE11690T	CE11890C
LOR Indicator Syringe	✓	✓	✓	✓	✓	✓	✓	✓
Epidural Needle 16G, 90mm	✓	✓	✓		✓	✓	✓	
Epidural Needle 18G, 90mm				✓				✓
Spinal Needle 25G, 123mm	✓	✓	✓					
Spinal Needle 27G, 123mm				✓				
Epi. Catheter 19G	✓				✓			
Epi. Catheter- Reinforced 19G		✓				✓		
Epi. Catheter- Soft tip 19G			✓				✓	
Epi. Catheter 21G				✓				✓
Medicine Filter 0.22μ	✓	✓	✓	✓	✓	✓	✓	✓
5ml/10ml syringe	✓	✓	✓	✓	✓	✓	✓	✓

Notice:
80mm epidural needle is optional. Subcutaneous needle (18G,22G,25G) is optional.